





VECTASTAIN[®] ABC-AmP Reagent (Standard, Western Blot Detection)

AK-6000

Product Images



Short Description

The VECTASTAIN ABC-AmP is an amplified ABC alkaline phosphatase reagent for the detection of mouse or rabbit primary antibodies on nitrocellulose or PVDF membranes. The signal can be visualized using either a chemiluminescent / chemifluorescent or a chromogenic substrate. Used with the DuoLuX Chemiluminescent / Chemifluorescent Substrate, the VECTASTAIN ABC- AmP system produces a very high and sustained light emission signal, with low background, and permanent fluorescence.

Features:

- Uses biotin/streptavidin technology
- Can detect as little as 1 pg target protein when used with DuoLuX Substrate
- Standard Kit contains Reagent A and Reagent B only
- Sufficient reagents to develop twenty 100 cm² blots

Additional Information

Unit Size	1 kit
Applications	Blotting Applications, Elispot
Technology	ABC (Streptavidin/Biotin) Systems
Detection Enzyme(s)	Alkaline Phosphatase
Target Species	No antibody included
Format	Concentrate

