




WestVision™ Block and Diluent for Western Blots

SP-7000-500

Product Images



Short Description

WestVision™ Block and Diluent is a ready-to-use Tris buffer-based blocking solution intended for Western or dot blot applications. This solution is a proprietary formulation containing Tween 20 that improves signal intensity and resolution without increasing background. To achieve the best results, use WestVision™ Block and Diluent for all blocking and diluting steps including primary antibody and secondary antibody detection reagents. WestVision™ Block and Diluent is compatible with both alkaline phosphatase- and peroxidase-based detections systems, as well as with both chemiluminescence and chromogenic development.

Additional Information

Unit Size	500 ml
Applications	Blotting Applications, Elispot
Blocking Action	Non-Specific Protein Blocking
Format	Ready-to-Use

