



# Fluorescein Rabbit Anti-Goat IgG Antibody

## Product Images



## Short Description

---

### Features:

- Affinity-purified, ultrapure, high affinity antibody
- Suitable for use with goat, sheep, and bovine IgG primary antibodies
- Thoroughly adsorbed against serum and immunoglobulins from potentially interfering species
- Unless otherwise specified, antibodies recognizes both heavy and light chains (H+L)
- Optimally labeled with fluorescein to provide the brightest label for fluorescence microscopy
- Supplied in solution
- Excitation: 495 nm
- Emission: 515 nm
- Color: Green

### Additional Information

---

Unit Size	1.5 mg
Applications	Immunofluorescence, In situ hybridization, Blotting Applications, Flow Cytometry/Cell Separation
Target Species	Goat
Conjugate	Fluorescein
Color of Fluorescence	Green
Host Species	Rabbit
Format	Concentrate

