



# Normal Horse Serum Blocking Solution, 2.5%

## Product Images



## Short Description

---

2.5% Normal Horse Serum (S-2012) is a ready-to-use blocking solution that is prepared identically to our undiluted normal serum. It can be used for blocking non-specific binding or as an antibody diluent.

### Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Same serum preparation that is included in the ImmPRESS® Detection Systems and the R.T.U. VECTASTAIN® ABC Kits

## Additional Information

---

Unit Size	50 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence
Blocking Action	Non-Specific Protein Blocking
Format	Ready-to-Use

