



# ImmPRESS® HRP Universal (Horse Anti-Mouse/Rabbit IgG) PLUS Polymer Kit, Peroxidase

## Product Images



## Short Description

---

The ImmPRESS polymerized reporter enzyme staining system uses novel conjugation and micropolymer chemistries to create a highly sensitive, ready-to-use, one-step detection system. The inclusion of BLOXALL® endogenous enzyme quenching solution and ImmPACT® DAB EqV substrate in the PLUS kit provides added convenience with optimized, matched reagents in one product.

### Features:

- High sensitivity and low background
- Enhanced accessibility to nuclear and membrane antigens
- Ready-to-use in a convenient dropper bottle
- Shorter assay time
- Simplified multiple labeling
- Designed for use with both rabbit and mouse primary antibodies

The ImmPRESS (15 ml) Kit will stain approximately 75-150 sections.

**Note:** This polymer detection kit should not be used to stain tissues from rodents or rabbits due to its reactivity with IgG present in the tissues of these species.

## Additional Information

---

Unit Size	1 Kit
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization, ELISAs
Technology	Micropolymer Reagents
Detection Enzyme(s)	Peroxidase
Conjugate	Micropolymer HRP
Target Species	Universal (Mouse/Rabbit)
Format	Ready-to-Use

