



NanoLink[®] 4FB Magnetic Beads 1.0 μm


M-1001

Product Images

Additional Information

Bead Size	1 μm
Loading Capacity	≥ 25 nmol/mg
Reactivity	HyNic
Applications	Antibody Labeling, Aptamers, In vitro Diagnostic (IVD), Laboratory Developed Tests (LDT), Analyte Specific Reagents (ASR), Next-Generation Sequencing (NGS), Photocrosslinking Studies
Recommended Storage	4°C

Products in this set

	NanoLink[®] 4FB Magnetic Beads 1.0 μm, 1.00 mL
--	--

