



# MagnaLink<sup>®</sup> Amino Magnetic Beads 2.8 µm

## M-1005

[Product Images](#)

### Additional Information

Bead Size	2.8 µm
Loading Capacity	≥ 25 nmol/mg
Reactivity	NHS Ester
Applications	Antibody Labeling, Aptamers, In vitro Diagnostic (IVD), Laboratory Developed Tests (LDT), Analyte Specific Reagents (ASR), Next-Generation Sequencing (NGS), Photocrosslinking Studies
Recommended Storage	4°C

### Products in this set



MagnaLink<sup>®</sup> Amino Magnetic Beads 2.8 µm, 1.00 mL

