



# Sulfo S-4FB Linker (Water Soluble)

## S-1008

### Product Images



## Short Description

---

The S-4FB heterobifunctional crosslinker is fundamental to the SoluLINK bioconjugation technology. Sulfo-S-4FB reacts with primary amines on proteins (through lysine amino acids) or amino-modified oligonucleotides or surfaces, introducing a 4FB linker that forms stable covalent conjugates with biomolecules possessing HyNic linkers. Sulfo-S-4FB contains a Sulfo-NHS ester, which makes this linker water soluble. No DMF or other organic solvents are necessary to dissolve the linker prior to use.

## Additional Information

---

Label/Modifier Type	4FB
Reactivity	HyNic
Recommended Storage	Desiccated: -15° to -25°C
Applications	Antibody Labeling, Antisense/RNAi, Aptamers
Other Name(s)	SoluLINK Bioconjugation

## Products in this set

---



**Sulfo S-4FB Linker (Water Soluble)**



**Sulfo S-4FB Linker (Water Soluble)**

