



Anhydrous DMF

S-4001-005

[Product Images](#)



Short Description

Anhydrous DMF (Dimethylformamide) is used for solubilizing S-HyNic, S-4FB, and other SoluLINK linkers. Anhydrous DMF is packaged with molecular sieves to minimize water contamination and preserve quality.

Additional Information

Recommended Storage	15° - 25°C
Applications	In Situ Proximity Ligation
Other Name(s)	Dimethylformamide SoluLINK Bioconjugation

Products in this set



Anhydrous DMF, 5 x 1.50 mL

