



Conjugation Buffer (10X)

S-4002-005

[Product Images](#)



Short Description

Conjugation Buffer is a high buffering capacity solution that is optimized for the conjugation reaction between HyNic and 4FB groups. Proteins, oligonucleotides, and other biomolecules are buffer exchanged into Conjugation Buffer after modification with S-HyNic and S-4FB, and prior to conjugation, to ensure efficient conjugation of these linkers.

Additional Information

Recommended Storage	2° - 8°C – Do Not Freeze
Applications	In Situ Proximity Ligation
Other Name(s)	SoluLINK Bioconjugation

Products in this set



Conjugation Buffer (10X), 5 x 1.50 mL

