



# Zeba™ Desalting Columns (0.5 mL, 7kDa)

## S-4024-010

[Product Images](#)

**SoluLINK®**

bioconjugation products  
and services are now  
only available from  
Vector Laboratories



## Short Description

---

Zeba™ Desalting Columns contain a high performance resin that offers exceptional desalting and protein recovery characteristics compared to other commercially available resins. These 0.5 mL Zeba Columns each have a capacity of up to 130 µL of material. Samples containing as little as 0.25 mg/mL of protein can be processed, providing exceptional protein recovery and > 95% retention of salts and other small molecules < 7 kDa. Zeba Desalting Columns are recommended for processing proteins and other biomolecules with a molecular weight of > 25 kDa to ensure maximum recovery.

Zeba is a trademark of Thermo Fisher Scientific.

## Additional Information

---

Recommended Storage	2° - 8°C – Do Not Freeze
Applications	In Situ Proximity Ligation
Other Name(s)	SoluLINK Bioconjugation

## Products in this set

---



**Zeba Desalting COLUMNS, 10.00 each**

