



# ImmEdge<sup>®</sup> Hydrophobic Barrier PAP Pen

## Product Images



## Short Description

---

The ImmEdge Pen (H-4000) is an improvement over the PAP pen for immunohistochemistry, immunocytochemistry, immunofluorescence, and *in situ* hybridization. It is designed to provide a heat-stable, water-repellent barrier that keeps reagents localized on tissue specimens and prevents mixing of reagents when multiple sections are mounted on the same slide.

### Features:

- Light blue barrier
- Compatible with enzyme or fluorescence-based detection systems
- Soluble in all commonly used clearing agents

## Additional Information

---

Unit Size	2 pens
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization
Safety Title	Cancer and Reproductive Harm

