



10X CASEIN SOLUTION

SKU: SP-5020-250



DESCRIPTION

10x Casein Solution is a general blocking agent for histochemical, nucleic acid and protein blotting, and other applications. It is supplied as 250 ml of a 10x concentrate.

SPECIFICATIONS

Blocking Action	Non-Specific Protein Blocking
Format	Concentrate
Unit Size	250 ml
Storage	2-8 °C
Instructions	

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Usage Summary

This product is designed to be used as a general blocking agent/diluent for nucleic acid and protein blotting applications. Prior to use, 10x Casein Solution should be diluted to 1x with distilled or deionized water. For prolonged storage of diluted Casein Solution, sodium azide should be added to give a final concentration of 0.08% (w/v) and the diluted solution refrigerated.

Applications

Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.