



NORMAL GOAT SERUM BLOCKING SOLUTION, 2.5%

SKU: S-1012-50



DESCRIPTION

2.5% Normal Goat Serum can be used for blocking non-specific binding or as an antibody diluent.

Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Ready-to use 2.5% solution (v/v) in buffer containing 0.08% sodium azide
- Same serum preparation that is included in the ImmPRESS[®] Detection Systems

SPECIFICATIONS

Blocking Action Non-Specific Protein Blocking

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Format	Ready-to-Use
Unit Size	50 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.