



# **NORMAL GOAT SERUM BLOCKING SOLUTION, 2.5%**

**SKU:** S-1012-50



---

## **DESCRIPTION**

2.5% Normal Goat Serum can be used for blocking non-specific binding or as an antibody diluent.

### **Features:**

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Ready-to use 2.5% solution (v/v) in buffer containing 0.08% sodium azide
- Same serum preparation that is included in the ImmPRESS<sup>®</sup> Detection Systems

## **SPECIFICATIONS**

**Blocking Action** Non-Specific Protein Blocking

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



<b>Format</b>	Ready-to-Use
<b>Unit Size</b>	50 ml
<b>Applications</b>	Immunohistochemistry / Immunocytochemistry, Immunofluorescence

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**