



## **NORMAL HORSE SERUM BLOCKING SOLUTION, 2.5%**

**SKU:** S-2012-50



---

### **DESCRIPTION**

2.5% Normal Horse Serum (S-2012) is a ready-to-use blocking solution that is prepared identically to our undiluted normal serum. It can be used for blocking non-specific binding or as an antibody diluent.

#### **Features:**

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Same serum preparation that is included in the ImmPRESS<sup>®</sup> Detection Systems and the R.T.U. VECTASTAIN<sup>®</sup> ABC Kits

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## SPECIFICATIONS

**Blocking Action** Non-Specific Protein Blocking

**Format** Ready-to-Use

**Unit Size** 50 ml

**Applications** Immunohistochemistry / Immunocytochemistry, Immunofluorescence

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**