



NORMAL HORSE SERUM BLOCKING SOLUTION, 2.5%

SKU: S-2012-50



DESCRIPTION

2.5% Normal Horse Serum (S-2012) is a ready-to-use blocking solution that is prepared identically to our undiluted normal serum. It can be used for blocking non-specific binding or as an antibody diluent.

Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Same serum preparation that is included in the ImmPRESS® Detection Systems and the R.T.U. VECTASTAIN® ABC Kits

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

Blocking Action Non-Specific Protein Blocking

Format Ready-to-Use

Unit Size 50 ml

Applications Immunohistochemistry / Immunocytochemistry, Immunofluorescence

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.