



GOAT ANTI-RABBIT IGG ANTIBODY (H+L)

SKU: CI-1000-1.5



DESCRIPTION

AMCA Goat Anti-Rabbit IgG Antibody can be used for immunofluorescence and other applications. Optimal F/P ratios have been established for each conjugate to ensure maximum fluorescence with minimal background staining.

Features:

- Recognizes both heavy and light chains (H+L)
- Optimally labeled with AMCA to provide the brightest label for fluorescence microscopy
- Supplied in solution
- Excitation: 350 nm
- Emission: 450 nm
- Color: Blue

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

Color of Fluorescence	Blue
Format	Concentrate
Formulation	10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide
Maximum Emission	448-454 nm.
Maximum Excitation	345-355 nm
Unit Size	1.5 mg
Storage Instructions	2-8 °C
Usage Summary	The recommended concentration range for use is 5-20 µg/ml
Applications	Immunofluorescence, In situ hybridization, Flow Cytometry/Cell Separation
Target Species	Rabbit
Concentration	1.5 mg active conjugate/ml
Conjugate	AMCA
Reactive Species	Goat
Source Species	Rabbit
Host Species	Goat

TECHNICAL INFORMATION

The goat anti-rabbit Ig antibodies are prepared by hyperimmunizing animals in a manner that produces high affinity antibodies. These are then purified by an affinity chromatography procedure designed to remove any low affinity antibodies which may be present. Cross-reactivities that are likely to interfere with specific labeling are removed by solid-phase adsorption techniques. The final product is then subjected to rigorous quality control assays including immunodiffusion, solid-phase enzyme immunoassays, gel electrophoresis and solid-phase binding assays. In preparing the labeled antibodies, great care is taken to ensure the maximum degree of labeling with no alteration in the specificity and affinity of the antibody. The labeled antibody then undergoes a further series of quality control assays, including immunohistochemical analysis.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.