




AMINO-DPEG®₂₄-ACID

SKU: QBD-10317



DESCRIPTION

Amino-dPEG®₂₄-acid, product number QBD-10317, is a PEGylated amino acid. A primary amine group and a propionic acid group sit on opposite ends of a long (76 atoms), single molecular weight polyethylene glycol (PEG) spacer of discrete length (dPEG®). The highly hydrophilic, discrete PEG spacer (Đ = 1) increases the hydrodynamic volume and imparts water solubility to conjugates that incorporate it. Moreover, Amino-dPEG®₂₄-acid benefitted numerous applications, including peptide synthesis, peptide modification, surface modification, oligonucleotide sequencing, drug delivery, and small molecule modification.

SPECIFICATIONS

CAS Number	756526-07-4
Molecular Weight	1146.35; single compound
Chemical Formula	C ₅₁ H ₁₀₃ NO ₂₆
Purity	> 97%
Unit Size	100 mg, 1000 mg
Solubility	Methylene chloride, DMAC or DMSO or water.
Spacers	dPEG® Spacer is 76 atoms and 89.0 Å

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions


-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.

Shipping Instructions Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.