



ANHYDROUS DMF

SKU: S-4001-005



DESCRIPTION

Anhydrous DMF (Dimethylformamide) is used for solubilizing S-HyNic, S-4FB, and other SoluLINK linkers. Anhydrous DMF is packaged with molecular sieves to minimize water contamination and preserve quality.

SPECIFICATIONS

Unit Size	5 x 1.5 mL
Storage Instructions	15° - 25°C
Applications	In Situ Proximity Ligation

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



CITATIONS



Powered by Bioz © 2023 See more details on Bioz

LEGAL NOTICE

Products are for research use only, not for use in diagnostic or therapeutic procedures or for use in humans. Products are not for resale without express written permission of Seller. No license under any patent or other intellectual property right of Seller or its licensors is granted or implied by the purchase unless otherwise provided in writing.

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.