



ANIMAL-FREE BLOCKER, 5X CONCENTRATE

SKU: SP-5030-250



DESCRIPTION

Animal-Free Blocker (5x) is a plant-derived agent and diluent for immunohistochemistry, nucleic acid and protein blotting, and other applications. This reagent contains no animal-derived protein and can be used as an alternative to sera, BSA, casein, or non-fat dry milk. It is supplied as 250 ml of a 5x concentrate.

SPECIFICATIONS

Blocking Action	Non-Specific Protein Blocking
Unit Size	250 ml
Storage	2-8 °C
Instructions	

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Usage Summary

This product is designed to be used as a general blocking agent/diluent for nucleic acid and protein blot applications. This reagent contains no animal-derived proteins and can be used as an alternative to sera, BSA, casein or non-fat dry milk in situations when animal protein-based blocking solutions cause increased background or interfere with blot development. Prior to use, 5x Animal-Free Blocker™ should be diluted to 1x with distilled or deionized water. For prolonged storage of diluted Animal-Free Blocker™, sodium azide should be added to give a final concentration of 0.08% (w/v) and the diluted solution refrigerated.

Applications

Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications, Elispot

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.