



ANTI-DNP (DINITROPHENYL)

SKU: SP-0603-1



DESCRIPTION

Dinitrophenyl (DNP) is not found endogenously in tissues and is an excellent alternative to biotin for use in applications requiring biotin-free systems. Probes or antibodies labeled with DNP can be detected using alkaline phosphatase labeled DNP or using the biotinylated antibody followed by a streptavidin or avidin system.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide.
Unit Size	1 mg
Storage Instructions	2-8 °C; Store frozen for long term storage
Applications	In situ hybridization, Blotting Applications

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Concentration	1 mg active conjugate/ml
Conjugate	Unconjugated
Reactive Species	Rabbit
Host Species	Rabbit

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.