



GOAT ANTI-GSL I

SKU: AS-2104-1



DESCRIPTION

These highly purified specific antibodies provide excellent intermediate reagents for localizing lectin receptors in tissues and on cells or following lectin transport in neuronal tracing.

These antibodies to lectins are produced by immunizing goats with pure lectins. Following conventional purification steps, specific antibody is isolated by affinity chromatography on lectin-agarose columns. Each anti-lectin is supplied in solution.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
Unit Size	2 mg

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions 2-8 °C; Store frozen for long term storage

Usage Summary The recommended concentration range for use is 1 to 10 µg/ml.
When reconstituting, do not vortex.

Applications Immunohistochemistry / Immunocytochemistry,
Immunofluorescence, Blotting Applications, Elispot, ELISAs,
Glycobiology

Concentration 2 mg active conjugate/ml

Conjugate Unconjugated

Reactive Species Goat

Host Species Goat

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.