



RABBIT ANTI-PHASEOLUS VULGARIS AGGLUTININ (E+L), UNCONJUGATED

SKU: AS-2300-1



DESCRIPTION

These highly purified specific antibodies provide excellent intermediate reagents for localizing lectin receptors in tissues and on cells or following lectin transport in neuronal tracing.

These antibodies to lectins are produced by immunizing rabbit with pure lectins. Following conventional purification steps, specific antibody is isolated by affinity chromatography on lectin-agarose columns. This anti-lectin is supplied in solution.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Unit Size	1 mg
Storage Instructions	2-8 °C; Store frozen for long term storage
Usage Summary	The recommended concentration range for use is 1 to 10 µg/ml.
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications, Elispot, ELISAs, Glycobiology
Concentration	2 mg/ml
Conjugate	Unconjugated
Host Species	Rabbit

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.