



GOAT ANTI-RICINUS COMMUNIS AGGLUTININ (RCA) I & II

SKU: AS-2084-1



DESCRIPTION

These highly purified specific antibodies provide excellent intermediate reagents for localizing lectin receptors in tissues and on cells or following lectin transport in neuronal tracing.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
Unit Size	1 mg
Storage Instructions	2-8 °C; Store frozen for long term storage
Usage Summary	The recommended concentration range for use is 1 to 10 µg/ml.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications, Elispot, ELISAs, Glycobiology
Concentration	2 mg active conjugate/ml
Conjugate	Unconjugated
Host Species	Goat

TECHNICAL INFORMATION

Each anti-lectin is supplied in solution. These antibodies to lectins are produced by immunizing goats with pure lectins. Following conventional purification steps, specific antibody is isolated by affinity chromatography on lectin-agarose columns.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Goat Anti-Ricinus Communis Agglutinin \(RCA\) I & II](https://vectorlabs.com/products/anti-ricinus-communis-agglutinin-rca-i-ii/)

<https://vectorlabs.com/products/anti-ricinus-communis-agglutinin-rca-i-ii/>