



GOAT ANTI-WHEAT GERM AGGLUTININ (WGA), UNCONJUGATED

SKU: AS-2024-1



DESCRIPTION

These antibodies to lectins are produced by immunizing goats with pure lectins. Following conventional purification steps, specific antibody is isolated by affinity chromatography on lectin-agarose columns. This anti-lectin is supplied in solution. These highly purified specific antibodies provide excellent intermediate reagents for localizing lectin receptors in tissues and on cells or following lectin transport in neuronal tracing.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
Unit Size	1 mg

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions 2-8 °C; Store frozen for long term storage

Applications Elispot, Immunohistochemistry / Immunocytochemistry, Blotting
Applications, Immunofluorescence

Concentration 2 mg active conjugate/ml

Conjugate Unconjugated

Reactive Species Goat

Host Species Goat

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [User Guide](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.