



# GOAT ANTI-WHEAT GERM AGGLUTININ (WGA), UNCONJUGATED

**SKU:** AS-2024-1



---

## DESCRIPTION

These antibodies to lectins are produced by immunizing goats with pure lectins. Following conventional purification steps, specific antibody is isolated by affinity chromatography on lectin-agarose columns. This anti-lectin is supplied in solution. These highly purified specific antibodies provide excellent intermediate reagents for localizing lectin receptors in tissues and on cells or following lectin transport in neuronal tracing.

## SPECIFICATIONS

<b>Format</b>	Concentrate
<b>Formulation</b>	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
<b>Unit Size</b>	1 mg

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



**Storage Instructions** 2-8 °C; Store frozen for long term storage

**Applications**                   Elispot, Immunohistochemistry / Immunocytochemistry, Blotting Applications, Immunofluorescence

**Concentration**                 2 mg active conjugate/ml

**Conjugate**                       Unconjugated

**Reactive Species**           Goat

**Host Species**                 Goat

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [User Guide](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**