




ANTIGEN UNMASKING SOLUTION, TRIS-BASED

SKU: H-3301-250



DESCRIPTION

Vector Laboratories' Antigen Unmasking Solutions are highly effective at revealing antigens in formalin-fixed, paraffin embedded tissue sections when used in combination with a high temperature treatment procedure. This Tris-based solution (H-3301) is pH 9.0 and is supplied as a convenient 100x concentrated stocks sufficient for preparation of 25 L of working solution.

SPECIFICATIONS

Unit Size 250 ml

Applications Immunohistochemistry / Immunocytochemistry, Immunofluorescence

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



TECHNICAL INFORMATION

CITATIONS





Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.