



AVIDIN D

SKU: A-2000-10



DESCRIPTION

Avidin D is isolated from egg white by a specific affinity system which results in a homogeneous, ultrapure avidin with very low nonspecific binding characteristics. Avidin D is recommended for use in solid-phase assays and many histochemical applications when the technique of sequential addition of reagents ("sandwich" technique) is employed.

SPECIFICATIONS

Molecular Weight	68000
Unit Size	10 mg
Storage Instructions	2-8 °C
Usage Summary	Dissolve in 0.1 M sodium bicarbonate. The recommended concentration range for use is 5-50 µg/ml.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications

Immunohistochemistry / Immunocytochemistry,
Immunofluorescence, Elispot, ELISAs

Conjugate

Unconjugated

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.