




AZDYE 350 AZIDE

SKU: CCT-1267



DESCRIPTION

AZDye™ 350 is a blue-fluorescent, azide-activated probe that reacts with terminal alkynes via a copper-catalyzed click reaction (CuAAC). It also reacts with strained cyclooctyne via a copper-free click chemistry reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures.


SPECIFICATIONS

CAS Number	N/A
Molecular Weight	437.47
Appearance	Grey to yellow solid
Extinction Coefficient	19,000
Purity	>95% (HPLC)
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C. Desiccate
Spectrally Similar Dyes	Alexa Fluor® 350, CF® 350, DyLight 350, AMCA
Excitation/Emission Maximum	346/445 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




ABS/EM SPECTRA



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.