





AZDYE 532 MALEIMIDE

SKU: FP-1042



DESCRIPTION



AZDye™ 532 Maleimide is a water-soluble, photostable, bright yellow-fluorescent dye with pH insensitive emission from pH 4 to pH 10. The excitation of AZDye™ 532 is ideally suited for the frequency-doubled Nd:YAG laser line. AZDye™ 532 dye can be conjugated to a variety of

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



antibodies, peptides, proteins, tracers, and amplification substrates, and is often used for the generation of stable signals in imaging and flow cytometry.


Maleimide is the most popular sulfhydryl-reactive group for conjugating the dye to a thiol group on a protein, oligonucleotide thiophosphate, or low molecular weight ligand. The maleimide group specifically and efficiently reacts with reduced thiols (sulfhydryl groups, -SH) at pH 6.5 to 7.5 to form a stable thioether bond. The resulting conjugates exhibit brighter fluorescence and greater photostability than the conjugates of many other spectrally similar fluorophores.

SPECIFICATIONS

Extinction Coefficient	78,000 cm ⁻¹ M ⁻¹
Reactivity	Thiol (Cys)
Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C.
Spectrally Similar Dyes	Alexa Fluor® 532, DyLight® 532
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.