




AZDYE 555 MALEIMIDE

SKU: FP-1168



DESCRIPTION

488/532



Laser
line

TRITC



Common
filter set

555



Excitation
max

580



Emission
max

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



AZDye™ 555 Maleimide is a water-soluble, bright orange-fluorescent dye with excitation ideally suited for the 532 nm or 555 nm laser lines and visualized with TRITC (tetramethylrhodamine) filter sets. AZDye™ 555 conjugates of antibodies, peptides, and proteins are pH insensitive from pH 4 to pH 10. Proteins can be labeled with AZDye™ 555 dye at high molar ratios without significant self-quenching, leading to brighter conjugates and more sensitive detection. AZDye 555 conjugates ideally suited for detection of low abundance targets.


Maleimide is the most popular sulfhydryl-reactive reactive group for conjugating the dye to a thiol group on a protein, oligonucleotide thiophosphate, or low molecular weight ligand. The maleimide group specifically and efficiently reacts with reduced thiols (sulfhydryl groups, -SH) at pH 6.5 to 7.5 to form stable thioether bond. The resulting conjugates exhibit brighter fluorescence and greater photostability than the conjugates of many other spectrally similar fluorophores.

SPECIFICATIONS

Molecular Weight	955.09 (protonated)
Extinction Coefficient	158,000 cm ⁻¹ M ⁻¹
Reactivity	Thiol (Cys)
Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C.
Spectrally Similar Dyes	Cy5, DyLight® 649
Excitation/Emission Maximum	555/572 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.