





AZDYE 568 CADAVERINE

SKU: FP-1083



DESCRIPTION



AZDye™ 568 Cadaverine is a carbonyl reactive building block used to modify carboxylic groups in the presence of activators (e.g. EDC, or DCC) or activated esters (e.g. NHS esters) through a stable amide bond. It also can be employed as a polar tracer and as a reactive dye for labeling

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



proteins via a carboxylic acid moiety.


AZDye™ 568 is a bright, and highly photostable, orange-fluorescent probe optimally excited by the 568 nm laser line on the Ar-Kr mixed-gas laser. This probe is water-soluble and its fluorescence is pH independent over a wide pH range. The brightness and photostability of this dye are best suited to direct imaging of low-abundance targets.

SPECIFICATIONS

Molecular Weight	892.92 (TFA salt)
Extinction Coefficient	93,000 cm ⁻¹ M ⁻¹
Reactivity	Aldehyde, carboxylic acid, ketone
Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C.
Spectrally Similar Dyes	Alexa Fluor® 568, CF™ 568
Excitation/Emission Maximum	578/602 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.