




AZDYE 594 AZIDE

SKU: CCT-1295



DESCRIPTION

AZDye™ 594 Azide is a bright, red-fluorescent azide-activated probe that reacts with terminal alkynes via a copper-catalyzed click reaction (CuAAC). It also reacts with strained cyclooctyne via a copper-free click chemistry reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures. AZDye™ 594 is well suited for the 561 nm and 594 nm laser lines of the Ar-Kr mixed-gas laser.

AZDye™ 594 Azide dye structurally is identical to Alexa Fluor® 594 Azide. Its absorption/emission spectra is a perfect match to spectra of many other fluorescent dyes based on sulfonated rhodamine core, including CF®594 Dye, DyLight® 594 and Alexa Fluor® 594.


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

CAS Number	N/A
Molecular Weight	846.97 (protonated)
Appearance	Red solid
Extinction Coefficient	88,000
Purity	>95% (HPLC)
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C. Desiccate
Spectrally Similar Dyes	Alexa Fluor® 594, CF® 594, DyLight® 594
Laser Line	561 nm or 594 nm
Excitation/Emission Maximum	590/617 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA




For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.