




AZDYE 594 BIOTIN

SKU: CCT-1396



DESCRIPTION

AZDye 594 Biotin (Alexa Fluor® 594 Biotin equivalent) can be used for detecting and quantifying biotin binding sites of avidin, streptavidin or neutravidin. This reagent overcomes major shortcomings of commonly used Biotin-4-fluorescein – poor solubility in aqueous media and pH dependent absorbance and emission. AZDye 594 Biotin is a water soluble reagent and its red fluorescence is pH independent from pH 4 to pH 10. A flexible PEG3 spacer between biotin moiety and fluorescent tag minimize steric hindrance involved in binding to avidin, streptavidin or neutravidin.

SPECIFICATIONS


CAS Number	N/A
Molecular Weight	1123.32
Appearance	Blue solid
Extinction Coefficient	88,000
Purity	>95% (HPLC)
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, DMSO, DMF, MeOH
Storage Instructions	-20°C. Desiccate

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Spectrally Similar Dyes	Alexa Fluor® 594, Atto™ 594, CF™ 594 Dye, DyLight™ 594
Laser Line	532 or 594 nm
Excitation/Emission Maximum	590/617 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[AZDye 594 Biotin](#)

<https://vectorlabs.com/products/azdye-594-biotin/>