




AZDYE 594 TETRAZINE

SKU: CCT-1364



DESCRIPTION

A bright, red-fluorescent probe used for detection TCO-tagged biopolymers. AZDye 594 Tetrazine demonstrates exceptionally fast cycloaddition kinetics (up to $30\,000\text{ M}^{-1}\text{ s}^{-1}$) with *trans*-cyclooctenes (TCO) as the dienophile, the fastest kinetics ever reported for any bioorthogonal reaction. In applications such as in vivo cancer imaging or pre-targeted cell labeling studies where rapid reaction kinetics is a must AZDye 594 Tetrazine probe would of great value.


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

Molecular Weight	892.98
Appearance	Red solid
Extinction Coefficient	92,000
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, MeOH, DMSO, DMF
Storage Instructions	-20°C. Desiccate
Spectrally Similar Dyes	Alexa Fluor® 594, Atto™ 594, CF™ 594 Dye, DyLight™ 594
Excitation/Emission Maximum	590/617 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA



SELECTED REFERENCES

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




1. Marquard, A.N., *et al.* (2020). Expanding the Scope of Antibody Rebridging with New Pyridazinedione-TCO Constructs. *Bioconjug Chem.*, **31(6)**, 1616-23. [[PubMed](#)]

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.