




AZDYE 647 GLY-GLY-GLY

SKU: CCT-1556



DESCRIPTION

Water-soluble, substrate for sortase mediated labeling of proteins. Sortase catalyzes a transpeptidase reaction between a specific internal sequence of a protein and an amine group present on the N-terminus of triglycine recently has become an area of great interest. This method of labeling proteins has been denoted as “Sortagging”.

AZDye™ 647 Gly-Gly-Gly is a water-soluble, pH-insensitive from pH 4 to pH 10, far-red-fluorescent probe. Its excitation ideally suited for the 633 nm or 647 nm laser lines. AZDye™ 647 is structurally similar to Alexa Fluor® 647, and spectrally is almost identical to Cy5 Dye, Alexa Fluor® 647, CF® 647 Dye, or any other Cyanine5 based fluorescent dyes.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	1232.46
Appearance	Dark blue solid
Chemical Formula	C52H77N7O19S4
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C. Desiccate


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Shipping Conditions
Shipping Instructions

Ambient temperature
Ambient temperature


ABS/EM SPECTRA



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.