




AZDYE 680 GLY-GLY-GLY

SKU: CCT-1578



DESCRIPTION

Water-soluble, substrate for sortase mediated labeling of proteins. Sortase catalyzes a transpeptidase reaction between a specific internal sequence of a protein and an amine group present on the N-terminus of triglycine recently has become an area of great interest. This method of labeling proteins has been denoted as “Sortagging”.

AZDye 680 is a bright and photostable near-IR dye that is spectrally identical to Alexa Fluor® 680, IRDye® 680RD and Cy5.5. AZDye 680 is often used for small animal imaging applications in the 700 nm wavelength and provides excellent clearance profiles in animals. Flow cytometry is another common application for this dye. AZDye 680 ideally suited for the 633 nm laser line.

SPECIFICATIONS

Molecular Weight	1434.67
Molecular Formula	C ₆₃ H ₁₀₅ BrN ₁₀ O ₁₆ S ₃
Appearance	Dark green solid
Chemical Formula	C ₆₃ H ₁₀₅ BrN ₁₀ O ₁₆ S ₃
Extinction Coefficient	185,000
Unit Size	1 mg, 5 mg, 25 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C. Desiccate
Spectrally Similar Dyes	Alexa Fluor® 680, IRDye® 680RD, CF™ 680 Dye

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Excitation/Emission Maximum 678/701 nm


Shipping Conditions

Ambient temperature

Shipping Instructions

Ambient temperature


ABS/EM SPECTRA



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.