




AZIDE AGAROSE

SKU: CCT-1038



DESCRIPTION

Azide agarose is a 6% crosslinked agarose resin that is activated with azide functional groups for covalent capturing alkyne-tagged biomolecules. This product is adequate to work in batch or column purification (low pressure).

SPECIFICATIONS

Activation Level	≥ 10 μmol Azide groups per mL resin
Appearance	Off-white slurry
Bead Size	Spherical, 50-150 μm
Enrichment Target	Alkyne-tagged proteins
Format	Batch or column purifications
Formulation	50% slurry
Preservative	20% ethanol
Unit Size	2 ml, 25 ml
Storage Instructions	2-8C
Support	6% Cross-linked agarose
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.