




AZIDO-PEG3-MALEIMIDE

SKU: CCT-AZ107



DESCRIPTION

Azido-PEG3-Maleimide enables simple and efficient incorporation of azido moiety onto antibodies, cysteine-containing peptides and other thiol-containing molecules.

Azido-PEG3-Maleimide degrades quickly (hours) at room temperature. This product is provided as a two-component easy-to-use kit. The stock solution of Azido-PEG3-Maleimide is prepared with a straightforward mix-and-stir protocol in 30 min.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	369.37
Appearance	Colorless to slightly yellow oil
Chemical Formula	C ₁₅ H ₂₃ N ₅ O ₆
Unit Size	25 mg, 100 mg, 1000 mg
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




SELECTED REFERENCES

1. Allo, B., *et al.* (2018). Clickable and High-Sensitivity Metal-Containing Tags for Mass Cytometry. *Bioconjugate Chem.* , **29(6)**, 2028-38. [[PubMed](#)]

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [User Guide](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.