




AZIDO-PEG4-ACID

SKU: CCT-AZ102



DESCRIPTION

Azido-PEG4-Acid is a polyethylene glycol-containing reagent with a terminal carboxylic acid that can be reacted with primary amine groups in the presence of activators (e.g. EDC, or DCC) forming a stable amide bond. The hydrophilic PEG4 spacer increases solubility in aqueous media.

SPECIFICATIONS

CAS Number	1257063-35-6
Molecular Weight	291.3
Appearance	Colorless to slightly yellow oil
Chemical Formula	C11H21N3O6
Purity	>95% (HPLC)
Unit Size	100 mg, 1000 mg, 5 g
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Azido-PEG4-Acid](https://vectorlabs.com/products/azido-peg4-acid/)

<https://vectorlabs.com/products/azido-peg4-acid/>