




BIOTIN (LONG ARM) NHS

SKU: SP-1210-50



DESCRIPTION

Biotin (Long Arm) NHS, water soluble, is a sulfonated derivative of our Biotin (Long Arm) NHS. Like Biotin (Long Arm) NHS, this reagent is useful for attaching biotin to proteins or other molecules containing a primary amino group. The spacer arm allows full accessibility of attached biotins by avidin derivatives or complexes. Unlike Biotin (Long Arm) NHS, this derivative is completely water soluble, and does not require DMF or DMSO to dissolve the activated biotin.

SPECIFICATIONS

Molecular Weight	556.6
Unit Size	50 mg
Species	Protein
Storage Instructions	Frozen (desiccated). Storage in solution is not recommended.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Tag/Group Incorporated Biotin

Reactive Group Primary Amino Group

Safety Title WARNING: NA –www.P65Warnings.ca.gov

Source Species Protein

CITATIONS




Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.