




BIOTIN-NHS ESTER

SKU: CCT-B102



DESCRIPTION

Biotin-NHS Ester is an amine-reactive building block used for modification of amine-containing molecule in organic media. NHS-activated biotin is water-insoluble reagent, thus requires an organic solvent such as DMSO or DMF to be dissolved. Short spacer arm adds minimal mass to modified molecules.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	341.38
Appearance	White crystalline
Chemical Formula	C ₁₄ H ₁₉ N ₃ O ₅ S
Purity	>98%
Unit Size	100 mg, 1000 mg
Solubility	DMSO, DMF
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




Shipping Instructions

Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.