




BIOTIN-PEG3-AMINE

SKU: CCT-1187



DESCRIPTION

Biotin-PEG3-Amine is a carboxyl-reactive building block with enhanced solubility in aqueous media used for derivatizing carboxyl groups in the presence of activators (e.g. EDC, or DCC) or activated esters (e.g. NHS esters) through a stable amide bond.

SPECIFICATIONS

CAS Number	359860-27-8
Molecular Weight	418.55
Appearance	Off-white to slightly grey amorphous solid
Chemical Formula	C ₁₈ H ₃₄ N ₄ O ₅ S
Purity	>98%
Unit Size	100 mg, 500 mg, 1000 mg
Solubility	DMSO, DMF, MeOH
Storage Instructions	-20C
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.