


BIOTIN-PEG3-AMINE

SKU: CCT-1187



Description

Biotin-PEG3-Amine is a carboxyl-reactive building block with enhanced solubility in aqueous media used for derivatizing carboxyl groups in the presence of activators (e.g. EDC, or DCC) or activated esters (e.g. NHS esters) through a stable amide bond.

Specifications

Unit Size	100 mg, 500 mg, 1000 mg
Molecular weight	418.55
Chemical composition	C ₁₈ H ₃₄ N ₄ O ₅ S
CAS	359860-27-8
Solubility	DMSO, DMF, MeOH
Purity	>98%
Appearance	Off-white to slightly grey amorphous solid
Storage Conditions	-20C
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.