



BOVINE SERUM ALBUMIN (BSA)

SKU: B-2007-10



DESCRIPTION

This product is produced from our ultrapure immunohistochemical grade bovine serum albumin. This biotinylated protein can be used for measuring binding capacity of avidin or streptavidin matrices, as a reagent for amplifying biotin-avidin detection systems, or as a positive control in dot blot assays.

SPECIFICATIONS

Formulation	10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% azide.
Unit Size	10 mg
Storage Instructions	2-8 °C
Applications	Blotting Applications, Affinity Chromatography
Concentration	5 mg/ml

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Conjugate

Biotinylated

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.