




# **BOVINE SERUM ALBUMIN (BSA), IMMUNOHISTOCHEMICAL GRADE**

**SKU:** SP-5050-500



---

## **DESCRIPTION**

This ultrapure grade of bovine serum albumin (BSA) can be used as a diluent or a blocking agent. It is free of impurities present in other grades of BSA which can introduce artifacts or increase background in immunohistochemical staining, ELISAs, or blot development.

## **SPECIFICATIONS**

<b>Blocking Action</b>	Non-Specific Protein Blocking
<b>Format</b>	Lyophilized
<b>Unit Size</b>	500 mg
<b>Storage Instructions</b>	2-8 °C

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## Applications

Immunohistochemistry / Immunocytochemistry,  
Immunofluorescence, In situ hybridization, Blotting Applications,  
Elispot, ELISAs, Flow Cytometry/Cell Separation

## CITATIONS




Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**