




BOVINE SERUM ALBUMIN (BSA), IMMUNOHISTOCHEMICAL GRADE

SKU: SP-5050-500



DESCRIPTION

This ultrapure grade of bovine serum albumin (BSA) can be used as a diluent or a blocking agent. It is free of impurities present in other grades of BSA which can introduce artifacts or increase background in immunohistochemical staining, ELISAs, or blot development.

SPECIFICATIONS

Blocking Action	Non-Specific Protein Blocking
Format	Lyophilized
Unit Size	500 mg
Storage Instructions	2-8 °C

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications

Immunohistochemistry / Immunocytochemistry,
Immunofluorescence, In situ hybridization, Blotting Applications,
Elispot, ELISAs, Flow Cytometry/Cell Separation

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.