



CONJUGATION BUFFER (10X)

SKU: S-4002-005



DESCRIPTION

Conjugation Buffer is a high buffering capacity solution that is optimized for the conjugation reaction between HyNic and 4FB groups. Proteins, oligonucleotides, and other biomolecules are buffer exchanged into Conjugation Buffer after modification with S-HyNic and S-4FB, and prior to conjugation, to ensure efficient conjugation of these linkers.

SPECIFICATIONS

Unit Size	5 x 1.5 mL
Storage Instructions	2° - 8°C - Do Not Freeze
Applications	In Situ Proximity Ligation

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



LEGAL NOTICE

Products are for research use only, not for use in diagnostic or therapeutic procedures or for use in humans. Products are not for resale without express written permission of Seller. No license under any patent or other intellectual property right of Seller or its licensors is granted or implied by the purchase unless otherwise provided in writing.

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.