


CY5.5 AZIDE

SKU: CCT-1059




Description

The NIR-fluorescent Cy5.5 azide can be reacted with terminal alkynes via a copper-catalyzed click reaction (CuAAC). It also reacts with strained cyclooctyne via a copper-free “click chemistry” reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures.

Cy5.5 Azide is a bright and photostable near-IR probe that spectrally similar to Alexa Fluor® 680, DyLight® 680, and IRDye® 680 dye. The Cy5.5 Azide is water-soluble, hydrophilic dye often a reagent of choice for assay where minimal non-specific binding and exceptional brightness is required. The fluorescence of Cy5.5 Azide is pH insensitive from pH 4 to pH 10 and produces minimal autofluorescence of biological specimens in this region of the spectrum. Fluorescence of this long-wavelength Cyanine dye is not visible to the human eye but is readily detected by most imaging systems.

Abs/Em Spectra

For research use only. Not intended for animal or human therapeutic or diagnostic use.



Specifications

Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Abs/Em Maxima	678/694 nm
Extinction Coefficient	190,000
Spectrally Similar Dyes	Alexa Fluor® 680, IRDye® 680RD, DyLight® 680
Molecular weight	985.13
CAS	N/A
Solubility	Water, DMSO, DMF
Purity	>95% (HPLC)
Appearance	Blue solid
Storage Conditions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.